

# Rocket® FBS Kits with Sampling Wand



At Rocket Medical we know our customers are our greatest source of new ideas and product improvements.

Our latest FBS Kits incorporate your user feedback and will mean better visualisation, easier scalp puncture and sample collection all combined in a smaller tray to save on packaging and storage space.

## CONTENTS:

- Rocket Amnioscope and Amnilume™ Single Use Light Source,
- Rocket Sampling Wand - combined blade and capillary tube holder
- 5 x 10cm x-ray detectable gauze swabs - green.
- 6 swab sticks
- 2.5gms petroleum jelly
- sponge holding forceps
- high heparin capillary tubes
- spare double-ended capillary tube holder
- tube covers
- Packed in cartons of 10 units, for single use

- ◆ **RELIABLE, COMPATIBLE ANALYSIS:** Rocket FBS Kits are specifically designed for use with Radiometer® ABL™, Siemens® Rapidlab, Instrumentation Laboratory® GEM™, ROCHE OMNI, cobas, Nova Biomedical lactate and other analysers
- ◆ **MAXIMUM ILLUMINATION NOW COMES AS STANDARD:**  
To give you brighter, more reliable visualisation every kit contains the latest **Rocket Amnilume™ Light Source.**
- ◆ **EASIER SAMPLE COLLECTION:** The Rocket Sampling Wand combines blade and capillary tube holders to help solve the problem of relocating the puncture site when transferring between blade and tube holder. It is as easy to use a ball point pen.

Compatible Analyser	Recommended sample volume	Rocket® FBS Kit with Sampling Wand	NHS Order codes
<b>ABL™ 5, 500, 70, 77, 80 FLEX</b>	65-85µl	<b>R57024-00-SW</b>	<b>FFY056</b>
<b>ABL™ 700 Series</b> with tube volume selected	35µl	<b>R57018-00-SW</b>	<b>FFY055</b>
<b>ABL™ 800 FLEX</b>	>85µl	<b>R57012-00-SW</b> <b>R57013-00-SW</b> <b>R57024-00-SW</b>	<b>FFY053</b> <b>FFY054</b> <b>FFY056</b>
<b>ABL™ 90 FLEX</b>	65µl	<b>R57024-00-SW</b>	<b>FFY056</b>
<b>ABL™ 90 FLEX PLUS</b>	45µl	<b>R57018-00-SW</b>	<b>FFY055</b>
<b>ABL™ 800 FLEX</b> with FLEXMODE selected	35µl	<b>R57018-00-SW</b>	<b>FFY055</b>
<b>Werfen - Instrumentation Laboratory</b> <b>GEM™ 4000, 5000 Series</b>	65µl	<b>R57024-00-SW</b>	<b>FFY056</b>
<b>ROCHE® OMNI c, OMNI s</b> <b>cobas b 121, 123, 221 systems</b>	35µl 85µl	<b>R57018-00-SW</b> <b>R57024-00-SW</b>	<b>FFY055</b> <b>FFY056</b>
<b>SIEMENS® (Bayer) Rapidlab™ Series</b>	100µl	<b>R57026-00-SW</b>	<b>FFY057</b>
<b>NOVA Biomedical StatStrip Lactate</b>	6µl	<b>R57028-00-SW</b>	<b>N/A</b>